
JuTrack

Release 2.0

mehran Sahandi far

May 29, 2020

CONTENTS


1	JuTrack Social	3
1.1	Requirements	3
1.2	How to Install	3
1.3	How to Activate	4
1.4	Main Menu	15
1.5	Passive Monitoring:	19
1.6	Active Monitoring	20
1.7	Ecological Momentary Assessment	20
2	JuSERVER	21
2.1	API	21
2.2	Dashboard	21
3	Developers	23
3.1	SAVE TO DATABASE	23
3.2	Synchronization with server	23
3.3	foreground services	23
3.4	crashHandler	23
3.5	OnBoarding	23
3.6	server Interface	23
4	how to contribute:	25
5	Source Code:	27
6	How to Cite :	29
7	Contact Us:	31

JuTrack is a platform for remote assessments of an individual's data in situ using smartphones. JuTrack is more than a single application, you can create, manage a study and, share your collected data using one platform.

A wide range of measurements using motion sensors (e.g. acceleration, gyroscope) as well as social activities (e.g. time spent in social networks) are extracted from the device while it is being used. Recorded are activity types (e.g. walking, running), relative location, and ecological momentary assessments. you may perform **Active** and **Passive** data collection and adopt the open-source codes for further requirements.

JuTrack platform includes:

- Android application
- Python dashboard
- python server API



JUTRACK SOCIAL


1.1 Requirements

- **JuTrack** relay on **Google Play Service** . If your Phone does not support it (e.g. some models of **Huawei** phones), you may not be able to use all modules.
- You need to be in **Germany** play store, otherwise you may not be able to download the app. To change it [visit](#)

1.2 How to Install

To install the JuTrack Social application go to the Google Play Store and search for **JuTrack Social**,

Or click [here](#) to access JuTrack Social in Google Play Store.




download and Install the JuTrack Social, then open it, and follow the next steps:

1.3 How to Activate

First check your Internet connection, activation and joining a study require an Internet connection. Then click on the **An Einer Studie Teilnehmen** (Join a study) on the welcome page.




The QR-Code scanning page will be opened, to enable the scanning please accept camera permission. Now, you can scan the provided Qr-Code within 1 minute.




Note: JuTrack Social will ask you to enable the Internet connection if it is not.

Important: provided QR-codes are for **One Time** usage, Using for a second time or in case of any problem, you will be redirected to **Retry** page.


In the next step, accept the necessary permissions. First, *Location Service* permission will pop-up, please enable it by selecting **ALLOW**




Then, JuTrack Social will take you to the **Data Usage Access** page in the Android setting to request for permission, which is necessary for *Application Usage Statistics* module.



Just select JuTrack Social from a list and **enable** it.




Almost done!, Just come back to the JuTrack Social by pressing back. Now depending on the Android version or Model of device, JuTrack may request to put it in **White List**. This step will request your Operation System (Android) or your device battery optimization program to permit JuTrack Social to work on background without of need for optimization. This is a **vital** step to keep the JuTrack Social running smoothly.




Note: This step might be different for your device or version of Android.


for **Samsung** devices with **Android 8 or above**:




Open **unmounted apps** as follow:




Then select JuTrack Social from a list and click done!



For **Huawei** device: select the JuTrack Social from a list and **disable** it.



for **Xiaomi** devices :




select JuTrack from a list and **Enable** it.

Congratulation! you finished all the necessary steps to Install and activate the JuTrack Social Application. As you may realized, JuTrack Social will show a notification as long as its running on background!


1.4 Main Menu

To access the main menu of JuTrack click **three dots** on the right corner of main page.



1.4.1 Administrator Page

In this menu, you have access to some administrative option to control the behavior of JuTrack. select **Einstellungen** from main menu, then insert the admin password (to get this password please contact us).




Note: you can enable the Phone data for synchronization (may cost extra charge depending on your carrier or contract).

Warning: Most of the options will be assigned automatically by the server during activation. Recommended NOT TO CHANGE THEM.

1.4.2 How to Leave a Study

you can leave a current study by selecting **Studie verlassen** (leave a study), you need to be connected to the Internet. During the leaving process (after accepting it) all the collected data by JuTrack will be Synchronized with the Study server.



1.4.3 Manual Synchronization

By selecting this option all data collected by JuTrack will be Synchronized with Study servers.

Note: Only use when its necessary.

1.5 Passive Monitoring:

In passive monitoring modules, the assessment will continue without any interaction from users. The following passive monitoring modules will collect data over specified criteria or time period automatically.

1.5.1 Location Service

location service records the user's relative location in a specified time period (i.e. every 10 minutes) or whenever there is a change in the user's location. for an energy-efficient implementation JuTrack benefits [Google's location service API](#) also, to meet the requirements of security and privacy, the location transformed in 3D space randomly for each user. this transformation is hidden for users and study owners which can not be reverted.

1.5.2 User Activity Detection

It detects users mode of transportation, for an efficient implementation JuTrack benefits [Google's activity detection API](#) in this module, JuTrack is able to detect up to 7 different activities every specified time (i.e. every 5 minutes) or whenever a new activity detected. the list of activities is as follow:

- IN_VEHICLE
- ON_BICYCLE
- ON_FOOT
- STILL
- UNKNOWN
- TILTING
- WALKING
- RUNNING

Besides detected activities, we also save an accuracy of detected activity in percentage.

Note: we only record the detected activity with accuracy >75%

1.5.3 Application Usage Statistics

Application Usage Statistics keeps track of applications used by the user. It captures and records the amount of time each application was in the **foreground** since the **Previous Midnight**. Among different applications such as social media, texting, games and etc, it also tracks the time spent on **phone calls** or Short Message Service (SMS) text messaging. it records these statistics every specified time period (i.e. every 1 hour).

Important: this module **does not** have access or track any personal information such as the content of the application or phone number.

1.5.4 Sensors

supported sensors:

- Accelerometer
- Gyroscope
- Barometer
- Magnetometer
- Gravity
- Light
- Linear Accelerometer
- Rotation
- Proximity

1.6 Active Monitoring

will soon be released

1.7 Ecological Momentary Assessment

will soon be released

JUSERVER

2.1 API

2.2 Dashboard

DEVELOPERS

3.1 SAVE TO DATABASE

3.2 Synchronization with server

3.3 foreground services

3.4 crashHandler

3.5 OnBoarding

3.6 server Interface

HOW TO CONTRIBUTE:

SOURCE CODE:

HOW TO CITE :

CONTACT US:

Interested in using JuTrack in your study? contact us:

To contribute as a developer please check [Developers](#) page.